

1. Introduction

Inertial Measurement Unit (IMU) Devices

- The **Inertial Measurement Unit (IMU) Devices** are wearable devices used for measuring and collecting data from the body motions of patients
- IMU devices connect to the mobile application wirelessly using the Bluetooth connection




Figure1: MetaMotionR board and R case [1]




Figure2: MetaMotionR Sensor (MMR) [2]

Limitations of existing mobile applications

- Wearable monitoring systems would:
 - collect data from multiple IMU devices
 - connect to external sensors
 - visualize data
- **No existing application supports all three functionalities**

2. Objectives

- Collect data from multiple IMU devices and external sensors
- Implement algorithms to analyze data collected
- Visualize real-time data
- Make an app suitable for most types of wearable health monitoring

3. Method

The New App is implemented based on the Wearable Health Monitoring application **MetaBase** and incorporates application **MetaWear's** external sensors and graphs functionalities




Figure3: MetaWear logo [3]



Figure4: MetaBase logo [4]

4. Results

1. Home Page




Figure5: MetaWear Homepage [5]

Figure6: New App Homepage

2. External Sensors

- Added external sensors Gpio0 – Gpio5




Figure7: MetaWear DeviceConfig page [6]

Figure8: New App DeviceConfig page 1

Figure9: New App DeviceConfig page 2

3. Data Visualization

- Each sensor has a graph, and a graph can be enabled using the Graph checkbox on the DeviceConfig page
- General Graphs: graphs of Accelerometer, Accelerometer Magnitude, Gyroscope, and Magnetometer
- Gpio Graphs: graphs of external sensors (Gpio0 – Gpio5)




Figure10: MetaWear StreamMonitor page [7]





Figure11: New App StreamMonitor page 1 (General Graphs)




Figure12: New App StreamMonitor page 2 (Gpio Graphs)

5. Conclusion

- The new application incorporates the functionalities of existing applications and provides greater flexibility for collecting Wearable Health Monitoring devices' data and monitoring the health conditions of patients

References

- [1] : <https://mbientlab.com/store/rectangle-case/>
- [2] : MbientLab, Sphinx, & Read the Docs. (2021). MetaMotionR Sensor (MMR). MetaWear Tutorials. Retrieved August 6, 2022, from <https://mbientlab.com/tutorials/MetaMotionR.html>
- [3] : MbientLab. (2014). MetaWear (3.6.0) [Mobile application software]. Retrieved from https://play.google.com/store/apps/details?id=com.mbientlab.metawear.app&hl=en_US&gl=US
- [4], [5], [6], [7] : MbientLab. (2016). MetaBase (3.5.0) [Mobile application software]. Retrieved from https://play.google.com/store/apps/details?id=com.mbientlab.metawear.metabase&hl=en_US&gl=US